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breformance LiquidColdCuring

Keep out of children's reach. Avoid contact with eyes and skin. Do not inhale vapors.

- R-codes:
R 11 Highly flammable
R 20/21/22 Harmful by inhalation, in contact with skin and if swallowed
R 33 Danger of cumulative effects
R 37/38 Irritating to respiratory system and skin
R 43 May cause sensitization by skin contact
S-codes:
S 2 Keep out of children's reach
S 9 Keep container in a well-ventilated place
S 16 Keep away from sources of ignition – no smoking
S 36/37 Wear suitable protective clothing and gloves
S 60 This material and its container must be disposed of as hazardous waste

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Instructions for use

breformance LiquidColdCuring

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Indication range

breformance LiquidColdCuring is a mixing liquid for the use of breformance as a cold-curing resin. An amine-free and hence color-stable resin is obtained for the fabrication of temporary crowns and bridges.

2. Indication

Temporary crowns and bridges, long-term temporary restorations to be worn over a period of up to 6 months, extraoral (indirect) repairs, covering attachments in combined dentures, fixation of milled acrylic denture teeth on CoCr dentures.

3. Safety recommendations and hazard warnings

breformance LiquidColdCuring is highly flammable and contains methyl methacrylate, ethandiol-, 1,2 dimethacrylate N,N-dimethyl-p-toluidine.

3.1 Symbols

4. Storage and durability information

breformance LiquidColdCuring must be stored at a dark place at room temperature (23 °C). Protect against direct sunlight and heat. Observe the expiration date on the label.

5. Processing

5.1. Mixing ratio

Add 4 drops of liquid to 0.2 g of powder at a room temperature of 23° C. Mix breformance LiquidColdCuring to obtain a liquid consistency. Add a larger or small quantity of liquid into the powder depending on the room temperature.

Please note: breformance LiquidColdCuring Liquid contains methyl methacrylate. Irritates eyes, skin and respiratory organs. Product must not come into contact with skin to avoid possible sensitization. Work under an extraction system. breformance LiquidColdCuring is exclusively intended to be used for extraoral processing. After powder and liquid have been mixed evenly to obtain a creamy consistency, allow the material to swell approx. 1 minute; afterwards the material is ready for processing.

Please note: Bubbles or flowmarks in the material can be avoided by ensuring a correct mixing ratio and by avoiding the adherence of any material to the bottom of the mixing cup. The material container must be sealed after taking out some material.

Apply a thin coat of breformance Insulating Liquid/Isoplast ip 750 ml (REF 540 0101 9) onto the plaster model before mixing the self-curing (autopolymerizing) resin. Prior to fixing the silicone key (haptosil D REF 540 0118 0), apply breformance LiquidColdCuring Liquid to areas that are difficult to access to allow that the mixed material flows more easily into these areas.

The processing time span is approx. 3 to 5 minutes depending on the room temperature. The time for polymerizing at air is approx. 10 minutes. Faster and bubble-free curing is achieved by the use of a pressure polymerization pot (3 to 4 bar, water temperature of 40° C, 5 minutes). If the restoration is built up using several polymerization steps, grinding and wetting of the surface with breformance LiquidColdCuring (brush) is required before adding any additional material (porcelain).

5.2. Grinding and polishing

Grinding with tungsten carbide burs, (e.g. REF D194 KT 23, max. 15,000 rpm)

Prepolishing with Rodeo round brush, (REF 350 0096 0, Ø 18 mm, max. 15,000 rpm)

High luster polishing with cotton buff, (REF 350 0065 0, max. 15,000 rpm)

Please note: It is recommended to perform grinding and polishing work under an extraction system.

5.3. Fabrication of an individual restoration

Wax-up with biotec wax, grey, (REF 510 0061 0)

Silicone key produced with haptosil D, (REF 540 0118 0)

Light-curing die varnish, opaque, red, (REF 540 0010 4)

Mark the line of the enamel area on the dentine tooth mould and cut back using a tungsten carbide bur for the application of transparent/incisal material. Powders of ceramic stains mixed with enamel or transparent materials (porcelains) have proved to be useful to achieve special effects when individualizing restorations.

6. Appropriation of cervical and enamel materials to dentine materials

Cervical 1	A2/A3/B3/C3	Enamel 1	A2/A3/B3/C2
Cervical 2	A 3,5/D3	Enamel 2	A3,5/D3

7. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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